

GoDirect™ Implant Line Expands with 4.7mm Diameter

GoDirect™, the one-piece implant for totally edentulous cases, now includes a 4.7mm diameter for expanded versatility. The 4.7mmD connection is rounder where the sides intersect and designed for placement with the new GoDirect Insertion Tool (GDIT).


All GoDirect implants feature a platform compatible with Implant Direct's new GPS™ overdenture attachment system as well as Zest Anchor's LOCATOR®. Available in three body diameters, GoDirect also offers multiple selections to accommodate variations in soft tissue. While GoDirect is designed primarily for use with an overdenture attachment system, it may also be converted for use with unique screw-receiving and ball attachment components.

Effective immediately, the 4.7mmD GoDirect is available in three lengths (8, 10 or 13mm) with a 1.5mmL shoulder height and Soluble Blast Media (SBM) surface. The All-in-One Packaging includes a transfer and comfort cap along with the one-piece implant.

Incorporate this new GoDirect implant into your practice!

[Watch Tutorial](#)

[Shop Online Now](#)



All-in-One Packaging includes items for a one-stage surgical procedure and impression making