



## New Paralleling Tool Now Available

Our newly released paralleling tool (PAR) features enhanced performance with our narrowest step drills. The SD2.3, which has a 2.0mmD apical end and a final dimension of 2.3mmD, is recommended as the initial drill in the surgical protocol for each Implant Direct™ implant system. The SD2.8 follows the SD2.3 to gradually widen the osteotomy with its 2.3mmD apical end and 2.8mmD final dimension.

**Effective immediately**, the new paralleling tool is included in complete surgical trays and is available for individual purchase.

**Update your surgical tray today!**

**Shop Online Now**

